Steri-7® - Bringing simplicity and peace of mind to infection prevention and control.



S-7XTRA® - Trusted in Australia since 2014



Hospital Grade Biocidal Disinfectant Cleaner

Proudly made in Australia

Listed with Claims on the ARTG: 232011 and 232014

EFFECTIVE AGAINST BACTERIA, ENVELOPED AND NON-ENVELOPED VIRUSES, SPORES, MOULD, FUNGI, YEAST, TUBERCULOSIS, PARVOVIRUS AND COVID 19.

KIND TO SURFACES.

Clean and Kill - One product for <u>all</u> situations.

- No need for wet contact times starts to work in 30 seconds in dirty conditions and goes on killing.
- The Reactive Barrier Technology (RBT) which offers residual effectiveness on touched surfaces for 24 hours against bacteria, COVID-19, and yeast with a further 72-hour biostatic surface protection for bacteria and COVID-19.
- Near pH Neutral and food safe.
- No rinsing required. Free from heavy detergent residues leaving streak-free finish.
- Contains multi-actives to reduce the chance of pathogens building resistance.
- Non-Flammable, Non-corrosive and safe to use on any material or surface including fabrics, stainless steel, plastics, wood, glass, soft metals and medical equipment.

Contact our team who is waiting to assist you with samples, advice and connect you with your local distributor.

Our products have been stringently tested, again and again,

to ensure that they deliver the highest levels of efficacy and safety.

• Contact Kill Times - Tested in <u>Dirty Conditions</u>. Times are reduced in clean conditions.

Bacteria	Kill Time		Kill Time
Acinetobacter baumannii	1 min	Virus - Enveloped	
Bacillus cereus	5 mins	Avian Flu H5N1	5 mins
Bacillus subtilis	5 mins	Bovine viral diarrhea virus	5 mins
Burkholderia cenocepacia	1 min	Influenza A H1N1 (Swine Flu)	5 mins
Burkholderia gladioli	1 min	Hepatitis C	5 mins
Burkholderia multivorans	1 min	Herpes EHV1	1 min
Campylobacter jejuni	5 mins	Herpes EHV4	1 min
Corynebacterium	5 mins	Human Immunodeficiency Virus (HIV)	5 mins
Enterobacter aerogenes	5 mins	Influenza A H1N1 (Swine Flu)	5 mins
Enterococcus hirae	45 seconds	Middle East respiratory syndrome coronavirus (MERS-CoV)	5 mins
Enterococcus faecium	1 min	Norovirus	5 mins
Escherichia coli (E. coli)	30 seconds	SARS and Coronavirus Covid-19	5 mins
Klebsiella pneumoniae NDM-1	1 min		
Legionella pneumophilia	1 min	Yeast, Mould and Fungus	
Listeria monocytogene	30 seconds	Aspergillus brasiliensis (niger)	5 mins
Methicillin Resistant Staphyloccus aureus (MRSA)	1 min	Aspergillus fumigatus	5 mins
Neisseria	5 mins	Candida albicans	5 mins
Proteus vulgaris	1 min	Penicillium	5 mins
Psudomonas aeruginosa	45 seconds	Trycophyton sp (Atheletes foot, Tinea)	5 mins
Rhodococcus equi	1 min		
Salmonella typhimurium	30 seconds	Virus – Non-Enveloped	
Salmonella enteritidis	1 min	Parvovirus	30 mins
Salmonella choleraesuis	1 min	Poliovirus	5 mins
Salmonella enteritidis	1 min	Rotavirus	5 mins
Salmonella typhimurium	1 min		
Staphyloccus aureus (MRSA)	45 seconds	Tuberculosis - TB	
Steptococcus equi	1 min	Mycobacterium terrae (TB) – Tested at 50g soiling	5 mins
Stenotrophomonas maltophilia	1 min	Mycobacterium Avium complex	5 mins
Vancomycin-Resistant Enterococcus (VRE)	1 min		
		Spores	
		Clostridium difficile (C. Diff)	1 min
		Clostridium perfringens	5 mins

Directions for use:

Remove heavy soiling from surface.

Spray or squirt on and wipe with a clean damp cloth evenly. Leave to dry. Use S-7XTRA wipes from pack or tub, wipe surface evenly and leave to dry.

For concentrate, follow the instructions on the label for dilution or ask your representative for advice. Use as above once diluted. Do not mix with other cleaning products or chemicals.

WELCOME TO THE HOME OF BI SECURITY

Change the way you clean with

S-7XTRA®

One product – All situations

S-7XTRA® is available in the following formats:

CONCENTRATES		
<text></text>	S7XCON	 S-7XTRA Unscented Concentrate 5L. Dilutes to between 50L to 250L for High-Risk to Low-Risk settings. Low-foaming formula makes it ideal in cleaning and scrubbing machines, dispensers, and manual cleaning equipment. Long lasting protection without wet contact times. Providing flexibility and a low-cost solution. Available in single unit or carton of 4. Follow label instructions for easy dilution or ask your representative for advice.
<text><text><text></text></text></text>	S7XCONF	 S-7XTRA FRESH Concentrate 5L. Fresh Ginger & Green Tea Fragrance. Dilutes to between 50L to 250L for High-Risk to Low-Risk settings. Low-foaming formula makes it ideal in cleaning and scrubbing machines, dispensers, and manual cleaning equipment. Long lasting protection without wet contact times. Providing flexibility and a low-cost solution. Suitable for bathrooms, floors and general areas. Available in single unit or carton of 4. Follow label instructions for easy dilution or ask your representative for advice. Also available in 5L Ready to Use.

See options for Diluting Stations and dispensers.

Diluting and Dispensing

PRECISION

We work with Precision Flow to provide solutions to meet your needs.

The AccuPro line of proportioners delivers dilution accuracy regardless of fluctuations in water pressure. AccuPro pressure regulation technology eliminates dilution variance by controlling the amount of water allowed into the unit, never allowing it to "lean out" the desired dilution target.

re regulated

- Configurations for many applications- 2, 3, and 4 button, and select valve models available
- Better chemical performance–consistent correct dilution
- Time and labor savings—solutions are prepared accurately with reduced preparation time

Hydro Systems offers many ancillary products that can further enhance the performance of your

AccuPro[®] dispenser.

We can configure the perfect system for your most demanding applications.

AccuPro dispensers can be customized to meet your company's specific business strategies.

Contact us to discuss all options available.



AccuPro



Bio-Misting

What is Bio-Misting?

It is a method of dispensing chemicals and disinfectants via a fine mist which is mechanically sprayed from a bio-misting machine. These can be electric, or battery operated. They can also be hand-held or on trolleys.

Where has it been used?

S-7XTRA is used in Australia in all areas listed below currently with the exception of agriculture at this time.

Bio-Misting is widely used in healthcare, emergency services, childcare, aged care, sporting facilities, public and private transport, hotels and resorts, aquatic centres, retail centres and supermarkets, pest control, farming and poultry management, equine facilities, airports, correctional facilities and recently in heritage restoration works.

Models



AEMOOG - 5L Cordless - suitable for large spaces and outdoor areas.



AEBURE - 3L Electric - 6m cord. For indoor spaces.

Car Care



AEUPGUN - 12.8L Cordless Backpack - suitable for large spaces and outdoor areas.



AEBMHAND - 3L Rechargeable Handheld. Suitable for vehicles



S7X250CC - S-7XTRA Car Care Kit. 250ml spray with microfibre cloth.

WELCOME TO THE HOME OF BI SECURITY

Г

Change the way you clean with

S-7XTRA®

One product – All situations

READY TO USE SOLUTIONS		
<text></text>	S7XRTU5L	 S-7XTRA Unscented Ready to Use 5L. 5 Litre Ready to Use is formulated for High-Risk settings providing you with certainty. No diluting required. Available in 5L in single units or carton of 4. For use on all hard surfaces and safe on all materials and fabrics including specialist chairs and equipment.
<text><text><image/></text></text>	S7XFRTU5L	 S-7XTRA FRESH Ready to Use 5L Fresh Ginger & Green Tea Fragrance. 5 Litre Ready to Use is formulated for High-Risk settings providing you with certainty. No diluting required. Available in 5L in single units or carton of 4. For use on all hard surfaces and safe on all materials and fabrics including specialist chairs. Suitable for bathrooms, floors and general areas.
	S7XRTU750	 S-7XTRA Unscented Ready to Use 750ml Tigger. No diluting required. Formulated for High-Risk settings. Available in 750ml single units or carton of 12. Spray on and clean surface with clean cloth. Allow to dry naturally. No rinsing required.

	S7XRTU500 and S7XRTU500T	 S-7XTRA Unscented Ready to Use 500ml Squirt Bottle and 500ml Trigger Bottle – Unscented. No diluting required. Formulated for High-Risk settings. Available in 500ml single units or carton of 12. Squirt or spray on and clean the surface with clean cloth. Allow to dry naturally. No rinsing required.
WIPES		
<text></text>	S7XW200	 Formulated for High-Risk Settings providing certainty. Tub wipes are formatted to suit surfaces where low moisture is preferred. Available in Tub of 200 Wipes in single units or carton of 6. Wipe surface evenly cleaning as you go. Leave to dry naturally, no need for wet contact times. Suitable for touchscreens, keyboards, and moisture sensitive equipment.
	S7PW80	 Formulated for High-Risk Settings providing certainty. Pillow/flat Packs of 80 Wipes in single units or carton of 12. Pillow wipes offer a larger, wetter format for use on large surface areas. Wipe surface evenly cleaning as you go. Leave to dry naturally, no need for wet contact times. Suitable for stainless steel counters, examination beds, specialist chairs and all hard surfaces.

For use in the following areas:

All healthcare settings, emergency services, dental, aged care, childcare, hospitality and catering, transport, and accommodation.



MATERIAL COMPATIBILITY CHART

S-7XTRA Ready to Use formula contains the following actives:

- Benzalkonium Chloride 0.09% w/v
- Poly Hydrochloride 0.09% w/v
- Didecyldimethylammonium Chloride 0.147%

S-7XTRA is non-toxic, non-corrosive and safe to use on the following surfaces and materials

Material	Use
Fabrics – cotton, wool, nylon	
Vinyl	
Leather	
Synthetic Leather (Ambla)	Bed and chair coverings
Transfer coated polyurethane (PU) fabric	Wearable healthcare products,
(Dartex [®])	drapes, and barriers
Glass	
Touchscreens – including glass and plastic	
Ceramic and Tile	
Grout	
Wood – raw and varnished	
Laminate	
Painted surfaces	
Rubber	
Closed Cell foam	Air filters and specialty sponges
Plastics - Hard and Flexible	
Transparent plastics	
Polyvinyl fluoride (PVF) (Tedlar [®])	Aircraft construction
Polypropylene	
Polyurethane	
Polyethylene	
Polycarbonate	
PVC acrylic	
Fiberglass	
Gas permeable cellulose acetate butyrate (CAB)	Contact lenses
Metals	
Carbon Graphite	
Copper	
Aluminium	
Brass	
Titanium	
Nickel	
Chrome and Chrome plated	

S-7XTRA® Hand Hygiene from Steri-7® Proudly made in Australia

Hand Hygiene is the most important thing you can do to stop the spread of

infections. It is simple: Wash your hands as often as you can with soap and water.

When water is not close at hand; we recommend using **S-7XTRA® Hand Rub**.

Complies with W.H.O and EN standards for Hand Hygiene as antibacterial.

- Dermatologically Tested and non-perfumed.
- Alcohol Free and non-flammable.
- Contains natural moisturisers which protect skin cells from environmental stress, dehydration and that maintains skin integrity.
- Works in dirty conditions and remains effective for up to 3 hours.
- Safe to use in sensitive environments where *alcohol can be an issue*.

S-7XTRA® Hand Rub is available in the following formats:

FOR LEAR BORDERS BARDE	S7HR5L	 S-7XTRA Hand Rub 5L Refill Available in single units or carton of 4. Made to refill bottles or dispensers. Proudly Made in Australia.
	S7HR600	 S-7XTRA Hand Rub 600ml Available in single units or carton of 10. Desktop size. Proudly Made in Australia.
	S7HR75	 S-7XTRA Hand Rub75ml Available in single units or carton of 24. Personal travel size perfect for carry- on, handbag, or car. Stay safe at home, work, play and all the places in between. Proudly Made in Australia.



ABOUT STERI-7 PTY LTD

We are a 100% Australian owned company and have been operating since 2014.

Our solutions have been manufactured in Australia since 2018 and comply with all therapeutic goods administration regulations and testing criteria.

The S-7XTRA range of products have been used in many industries throughout Australia including:

- + Healthcare
- Emergency Services
- ✦ Aged Care
- ✤ Dental
- + Childcare
- + Education
- ✦ Transport
- Cleaning Services
- + Hospitality
- ✦ Mining
- ✦ Aviation
- Sports and Recreation
- Heritage Conservation
- Mould Remediation

We are proud to work directly and through our distributors who are experts in the fields they supply to and can assist in the transition to S-7XTRA through training, education, process and protocol development and implementation.

We are also proud to work with our partners who supply dispensing and dilution solutions to enhance and streamline processes within your facilities and our supply chain partners who service our needs within Australia and NZ.

Contact our team who is waiting to assist you with samples, advice and connect you with your local distributor.

WWW.STERI-7.COM.AU (02) 8209 3846 INFO@STERI-7.COM.AU



S-7XTRA® follows the Steri-7 Philosophy – **Break, Treat, Prevent Break** – The chain of infection by washing hands and **Treat** – Surfaces with S-7XTRA to

Prevent – The spread of infection to minimise risk.

S-7XTRA® offers a simpler, smarter and safer solution to your infection prevention protocols making them the smarter choice.