

S-7XTRA

Cleaner / disinfectant for non-invasive medical devices and environmental surfaces



S-7XTRA® is a combined cleaner/disinfectant that is suitable for multiple uses and is designed for exceptional materials compatibility* and the tri-effect of “Break, Treat, Prevent”.



MICROBIAL EFFICACY

S-7XTRA® is bactericidal, virucidal, fungicidal, yeasticidal, and sporicidal* on non-invasive medical devices and environmental surfaces



TESTING

S-7XTRA® has been rigorously tested for use in medical, institutional, food, and domestic environments using industry-standard test methods.



FORMULATION

Non-corrosive, non-viscous, and pH neutral, S-7XTRA® is designed to work in dirty conditions and contains triple active ingredients with different modes of action to reduce the risk of antimicrobial resistance.



REGULATORY COMPLIANCE

S-7XTRA® is included in the Australian Register of Therapeutic Goods as a *class IIb medical device disinfectant* (non-scented products) and as a *listed disinfectant* (non-scented and 'FRESH' products).

PRODUCT OPTIONS



FORMAT	SIZE	PRODUCT CODE
S-7XTRA® Concentrate	5L Container	S7XC0N
S-7XTRA® FRESH* Concentrate	5L Container	S7XFC0N
S-7XTRA® Ready-to-Use (RTU)	5L Container	S7XRTU5L
	750ml Trigger Spray	S7XRTU750
	500ml Squeeze Bottle	S7XRTU500
	500ml Trigger Spray	S7XRTU500T
S-7XTRA® FRESH* Ready-to-Use (RTU)	5L Container	S7XFRTU5L
S-7XTRA® Wipes	Pillow Pack x80	S7XPW80
	Pillow Pack x100	S7XPW100
	Tub x200	S7XW200

SIMPLER. SMARTER. SAFER.

Break the chain of infection with S-7XTRA

*Always refer to cleaning and disinfection instructions provided by the Medical Device Original Equipment Manufacturer before using S-7XTRA®.

*Against bacillus subtilis only on medical devices. S-7XTRA® is also mycobactericidal (including TB) and sporicidal (including against C. diff) on environmental surfaces. Refer to the 'Microbiocidal Activity' tables in the *Instructions For Use* for targeted pathogens and contact times.

*S-7XTRA® FRESH (ginger and green tea scented) is intended for use on environmental surfaces only.

Steri-7 Pty Ltd

🏠 20/8 Tilley Lane FRENCHS FOREST NSW 2086

☎ 02 8209 3846

✉ info@steri-7.com.au

🌐 www.steri-7.com.au